

DDPA2D4

Cobalt™ XT DR MRI SureScan™
dual chamber implantable
cardioverter defibrillator (ICD) –
technology: iATP; header: DF4



Cobalt™ XT ICDs are used to treat sudden cardiac arrest and abnormal heart rhythms.

Contact a representative

Description

Cobalt™ XT ICDs feature exclusive algorithms to manage AF, 40 J energy available on the first shock and every shock,¹ and extended longevity. These devices are enabled with BlueSync™ technology, allowing for tablet-based programming and app-based remote monitoring, and are compatible with an automated heart failure risk stratification tool.

1. Medtronic Cobalt™ XT HF Quad MRI SureScan™ Model DTPA2QQ device manual.

Not what you're looking for?

Find more information, including other available configurations, for [Cobalt XT™ MRI SureScan™ ICD](#).

Careers

Contact & Support

Order Products [↗](#)

 Medical Supplies [↗](#)

[Product Manuals](#) 

[Product Security](#) 

[Information About Proposition 65 for California Customers](#) 

[Digital Millennium Copyright Act \(DMCA\) Notification](#)

Medtronic
Engineering the extraordinary

[Privacy Notice](#)



[Terms of Use](#)

[Accessibility Statement](#)

Operational Headquarters 710 Medtronic Parkway,
Minneapolis, MN 55432-5604 USA

©2025 Medtronic

